

Neodent Helix GM Dental Implant System

Implant

Ø 3.5			Ø 3.75			Ø 4.0			Ø 4.3			Ø 5.0			Ø 6.0		
																	
	115€ Acqua	105€ NeoPoros		115€ Acqua	105€ NeoPoros		115€ Acqua	105€ NeoPoros		115€ Acqua	105€ NeoPoros		115€ Acqua	105€ NeoPoros		115€ Acqua	105€ NeoPoros
8.0	140.943	109.943	8.0	140.976	109.976	8.0	140.982	109.982	8.0	140.948	109.948	8.0	140.953	109.953	8.0	140.1009	109.1009
10.0	140.944	109.944	10.0	140.977	109.977	10.0	140.983	109.983	10.0	140.949	109.949	10.0	140.954	109.954	10.0	140.1010	109.1010
11.5	140.945	109.945	11.5	140.978	109.978	11.5	140.984	109.984	11.5	140.950	109.950	11.5	140.955	109.955	11.5	140.1011	109.1011
13.0	140.946	109.946	13.0	140.979	109.979	13.0	140.985	109.985	13.0	140.951	109.951	13.0	140.956	109.956	13.0	140.1012	109.1012
16.0	140.947	109.947	16.0	140.980	109.980	16.0	140.986	109.986	16.0	140.952	109.952	16.0	140.957	109.957			
18.0	140.988	109.988	18.0	140.981	109.981	18.0	140.987	109.987	18.0	140.989	109.989	18.0	140.990	109.990			
			20.0		109.1043	20.0		109.1046									
			22.5		109.1044	22.5		109.1047									
			25.0		109.1045	25.0		109.1048									

Cover Screw

			24€
GM Cover Screw	0 mm	2 mm	
	117.021	117.022	

Healing Abutment

GM Healing Abutments



31€

	Ø 3.3	Ø 4.5	Ø 5.5	Ø 6.5
0.8 mm	106.207	106.213		
1.5 mm	106.208	106.214	106.250	106.254
2.5 mm	106.209	106.215	106.251	106.255
3.5 mm	106.210	106.216	106.252	106.256
4.5 mm	106.211	106.217	106.253	106.257
5.5 mm	106.212	106.218	-	-

GM Customizable Healing Abutments (peek)



46€

	1.5 mm	2.5 mm	3.5 mm	4.5 mm	5.5 mm	6.5 mm
Ø 5.5	106.223	106.224	106.225	106.226	106.227	-
Ø 7.0	-	106.228	106.229	106.230	106.231	106.232

Neodent Helix GM PROSTHETICS

GM Titanium Base

GM Exact
Titanium
Base

48€

AH				4 mm		
GH	0.8 mm	1.5 mm	2.5 mm	3.5 mm	4.5 mm	
Ø 3.5	135.355	135.356	135.357	135.358	135.359	
Ø 4.5	135.367	135.368	135.369	135.370	135.371	

AH				6 mm		
GH	0.8 mm	1.5 mm	2.5 mm	3.5 mm	4.5 mm	
Ø 3.5	135.361	135.362	135.363	135.364	135.365	
Ø 4.5	135.373	135.374	135.375	135.376	135.377	

GM Titanium
Base for
Bridge

48€

AH						
GH	0.8 mm	1.5 mm	2.5 mm	3.5 mm	4.5 mm	
Ø 3.5	135.399	135.400	135.401	135.402	135.403	
Ø 4.5	135.404	135.405	135.406	135.407	135.408	



20 N.cm



20 N.cm



20 N.cm

GM Implant Exact Impression
Transfer Open Tray

Regular	Long
108.162	108.163

31€

48€

GM Temporary Abutment Crown

	0.8 mm	1.5 mm	2.5 mm	3.5 mm
Ø 3.5	118.344	118.345	118.346	118.347
Ø 4.5	118.348	118.349	118.350	118.351



20 N.cm

GM Implant
Exact
Scanbody

108.183

25€

GM Implant
Scanbody
Titanium N/A

108.207

49€

For lab



GM Analogue
Hybrid

Ø3.5/3.75	101.103
Ø4.0/4.3	101.089

24€



Implant
Scanbody

108.181

25€



Neo GM
Screw Long

116.292

21€

GM Multi-Unit Abutment / GM Exact Mini Conical Abutment

GM
Multi-unit
Abutment

72€

GH	0.8 mm	1.5 mm	2.5 mm	3.5 mm	4.5 mm	5.5 mm
0°	115.243	115.244	115.245	115.246	115.247	115.248
17°	-	115.275	115.276	115.277	-	-
30°	-	115.278	115.279	115.280	-	-
45°	-	115.281	115.282	-	-	-

For 0° use Multi-unit Driver | 32N.cm
For 17°, 30°, 45° use Neo Screwdriver | 20N.cm

32 N.cm



20 N.cm



10 N.cm



Multi-unit
Impression Transfer
Open Tray

108.176

31€

10 N.cm



Titanium
Coping

118.302

31€

10 N.cm



Hybrid
Coping

118.382

31€

10 N.cm



Multi-unit
Protection Cap

106.268

34€



Metal
Scanbody, Ti

108.207

49€

For lab



Analog

101.020

24€



Digianalog

101.092

24€



Scanbody

108.196

43€



GM Multi-unit
Coping Screw

116.269

11€

Screwdriver

36€

Neo
Screwdriver



Multi-unit

Contra-angle		Torque wrench	Manual	Driver
extra short	105.146	short	105.133	104.058
short	105.135	medium	105.132	104.060
long	105.160	long	105.157	104.070

Novaloc

116€

GM Novaloc

GH	0.8 mm	1.5 mm	2.5 mm	3.5 mm	4.5 mm	5.5 mm
0°	102.161	102.162	102.163	102.164	102.165	102.166
15°	102.167	102.168	102.169	102.170	102.171	102.172

20 N.cm



20 N.cm



Impression Cap
(4 pcs)

2010.722 - NOV

106€


For lab



Analog (4 pcs)

2010.721 - NOV

106€



Processing
Package titanium
(2 pcs)

2010.601 - NOV

72€